

CRYOALFA SUPER

Super and Super C pag. 5/5

Super

Art.no. CA-S

Liquid freezing therapy



SUPER device for "Liquid freezing" treatment for precise application in daily dermatological indications

5-8µ gold-plated devices and applicators

cartridges with in-built valve filled with sterile liquid- N₂O gas (nitrous oxide, cryogenic gas with constant temperature of -89 °C / -127 °F) medical quality, filtered to 5µ

the standard applicator is interchangeably with the contact applicators 3,5 + 7mm width

SUPER Contact device using for daily medical-aesthetically indications preferably treatment in the epidermis, as for pigmentations, spots, haemangiomas

all contact applicators are interchangeably with the standard liquid-applicator

Super C

Art.no. CA-S-C

Contact freezing therapy



protection cap
Art.no.CAS-PCa



Applicator for liquid freezing therapy with dosing feeder
Art.no.CAS-DST



connecting to cartridges with in-built valves



Art.no.CA-KV 8
Art.no.CA-KV 12
Art.no.CA-KV 16
Art.no.CA-KV 25

Packaging



Super and Super contact in blister pack



Super and Super contact in case packed with 4 cartridges



contact area of 3mm Ø
Art.no. CAS-C3



contact area of 5mm Ø
Art.no. CAS-C5



contact area of 7mm Ø
Art.no. CAS-C7